

# Patient Return Electrode Cord




**E0560** 1 EACH

**Patient Return Electrode Cord/Clamp, Reusable**

Compatible with with RECQM generators  
 For use with E7508, E7509 and E7509-B patient return electrodes

- REM™ pin connector
- 4,6 m (15 ft) cord



**E0560SI** 1 EACH

**Patient Return Electrode Cord/Clamp**

For use with all electrosurgical generators equipped with any monitoring system for dual patient return electrode (Short Isolation)

- 4,6 m (15 ft) cord



**E0560E** 1 EACH

**Patient Return Electrode Cord/Clamp**

For use with E7509 and E7509-B patient return electrodes and ERBE generators

- 4,6 m (15 ft) cord



**E0560NCS** 1 EACH

**Patient Return Electrode Cord/Clamp**

For use with all electrosurgical generators which are not equipped with any monitoring system for dual patient return electrode (Non-contact Surface)

- 4,6 m (15 ft) cord



**E0560LI** 1 EACH

**Patient Return Electrode Cord/Clamp**

For use with all electrosurgical generators equipped with any monitoring system for dual patient return electrode (Large Isolation)

- 4,6 m (15 ft) cord

# Electrosurgical Accessories, Non-Sterile, Reusable Bipolar Cords



**E360130** 1 EACH

**Foot Switching Forceps Cord**

For use with ERBE/STORZ/WOLF generators only

Compatible with E700xxx or 700xxx bipolar forceps

- Moulded plug
- 4,5 m (15 ft)



**E360150L, E360150<sup>2</sup>** 1 EACH

**Foot Switching Bipolar Cord**

For use with Valleylab™ generators only

Compatible with E700xxx or 700xxx bipolar forceps

- Moulded plug
- 4,5 m (15 ft)



**E360131L** 1 EACH

**Foot Switching Forceps Cord**

For use with ERBE/STORZ/WOLF generators only

Compatible with E700xxx or 700xxx bipolar forceps

- Moulded plug
- 4,5 m (15 ft)



**92067** 1 EACH

**Reger™ Cable for Ethicon PowerStar Bipolar Scissors for International 2-pole Models, 5,0 m**

- Generator side plug: international 2-pole model
- Instrument side plug: Ethicon PowerStar Bipolar Scissors
- Mechanical properties: to connect Ethicon PowerStar Bipolar Scissors with generators from international models
- Electrical properties: according to IEC 60601-1, IEC 60601-2-2
- Reprocessing: 300 times autoclavable at 134 °C; holding time 10 min
- Combination items: Ethicon PowerStar Bipolar Scissors
- Colour: grey
- Surface: matt



**E360139L** 1 EACH

**Foot Switching Cord**

For use with Martin/Berchthold/Aesculap generators only

Compatible with E700xxx or 700xxx bipolar forceps

- 4,5 m (15 ft)



**E360135<sup>1</sup>** 1 EACH

**Foot Switching Cord**

For use with US generators only

Compatible with E700xxx or 700xxx bipolar forceps

- 4,5 m (15 ft)

1. Not allowed in some European countries

2. Short connector

# Electrosurgical Accessories, Non-Sterile, Reusable Bipolar Cords

**Steam sterilizable**



**E0020V** 1 EACH

**Foot Switching Bipolar Cord**

For use with the ForceTriad™ and Valleylab™ FX8 energy platforms

Also compatible with all Force™ series and SurgiStat™ II generators

For use with Valleylab™ bipolar foot switching forceps

- Moulded plug
- 4,5 m (15 ft)



**E0021S** 1 EACH

**Foot Switching Bipolar Cord**

For use with the ForceTriad™ and Valleylab™ FX8 energy platforms

Also compatible with all Force™ series and SurgiStat™ II generators

Compatible with most Storz bipolar instruments

- Moulded plug
- 4,5 m (15 ft)



**E0022W** 1 EACH

**Foot Switching Bipolar Cord**

For use with the ForceTriad™ and Valleylab™ FX8 energy platforms

Also compatible with all Force™ series and SurgiStat™ II generators

Compatible with most Wolf bipolar forceps

- Moulded plug
- 4,5 m (15 ft)

**Bipolar resection cords**



**FT0021S** 1 EACH

**Storz™ Bipolar Resection Cord**

- To be used with Storz™ bipolar Resectoscopes
- To be used with the ForceTriad™ energy platform
- Reusable
- 4,6 m cord



**FT0022W** 1 EACH

**Wolf Bipolar Resection Cord**

- To be used with Wolf bipolar resectoscopes
- To be used with the ForceTriad™ energy platform
- Reusable
- 4,6 m cord


# Electrosurgical Accessories, Non-Sterile, Reusable, Monopolar Cords



**E2999** 1 EACH

### Foot Switching Monopolar Laparoscopic Cord

- Accepts 4 mm male post connector
- 3 m (10 ft)




**360186** 1 EACH

### Foot Switching Monopolar Cord

For use with Erbe ACC/ICC generators

- 4 mm female connector
- 3,5 m (11,5 ft)
- Use with Monopolar Forceps 700XXX



**360187** 1 EACH

### Foot Switching Monopolar Cord

For use with Valleylab

- Jack 8 mm connector
- 3,5 m (11,5 ft)
- Use with Monopolar Forceps 700XXX

# Electrosurgical Accessories, Non-Sterile, Reusable Adapters



**E0017** 1 EACH

**Universal Monopolar Adapter**

Compatible with Valleylab™ and most electrosurgical generators for foot switching activation

- Accepts standard monopolar accessory cords up to 6 mm



**E0504-2** 1 EACH

**Non-REM™ Patient Return Electrode Adapter**

See Patient Return Electrode Adapter Quick Reference Index on page 37 for a list of compatible generators

- Accepts standard non-REM™ Adult PolyHesive™ patient return electrode (E7506)



**E0502-1** 1 EACH

**Monopolar Adapter**

Compatible with Valleylab™ and most electrosurgical generators for foot switching activation

- Accepts standard monopolar accessory cord pin connectors up to 3,2 mm



**E0507-B** 1 EACH

**Multiple Return/S Cord Adapter**

For use with all Valleylab™ REM-equipped generators

- Accepts Valleylab™ REM PolyHesive™ II patient return electrodes
- Enables use of two patient return electrodes
- Enables use of an S cord

**E0005-3C** 1 EACH

**Monopolar Adapter**

- Accepts 4 mm standard active accessory



**368587** 1 EACH

**Monopolar Adapter**

For use with Martin/ Berchtold generators

- Accepts Valleylab™ 3 pin electrosurgical pencils



**368687** 1 EACH

**Monopolar Adapter**

For use with Erbe ACC/ICC generators

- Accepts Valleylab™ 3 pin electrosurgical pencils