


Sydney IVF Embryo Transfer Catheter Sets

Used for the atraumatic transfer of embryos into the uterine cavity.

- Precurved guiding catheter facilitates catheter insertion.
- Atraumatic bulb tip eases passage through cervix.
- Soft, flexible transfer catheter minimizes endometrial trauma.
- Microvol™ technology reduces the volume of medium required to transfer embryos.
- Cervical stop comes set at 4 cm with an additional stop at 5 cm.



Global Product Number	Order Number	Guide Catheter Fr	Guide Catheter Length cm	Transfer Catheter Fr	Transfer Catheter Length cm	Description
Catheter Sets						
Silicone Hub						
G20683	K-JETS-6019-SIVF	6.6	16.7	2.8	23	catheters packed assembled
Polycarbonate Hub						
G18740	K-JETS-7019-SIVF	6.6	17.3	2.8	24	catheters packed separately
Components						
G48344	K-JET-7002-SIVF	malleable obturator to fit K-JETS-7019-SIVF				