



# Endo-Scrub<sup>®</sup> 2 Lens Cleaning Sheaths

FOR THE INTEGRATED POWER CONSOLE  
(IPC<sup>®</sup> SYSTEM)



Maximum Visualization  
at the Surgical Site

# Endo-Scrub® 2 Lens Cleaning Sheaths

FOR THE INTEGRATED POWER CONSOLE (IPC® SYSTEM)



Our Endo-Scrub® 2 Lens Cleaning Sheaths help ENT surgeons maintain maximum visualization at the surgical site without removing the endoscope. The software-driven Endo-Scrub 2 sheath delivers a high-power burst of irrigation to flush debris from the lens tip. When used correctly, the sheath provides outstanding and consistent view maintenance, and only adds 0.7 mm to the endoscope's barrel diameter profile.

## Improved Software and Sheath Design

Our newest Endo-Scrub 2 software is fully integrated with and optimized for our Integrated Power Console (IPC® System). Multiple software improvements ensure that all fluid is evacuated from the sheath to keep the lens from fogging. We also redesigned the sheaths to better secure the scope in the sheath and help prevent fluid seepage.



## Key Features


- Automatically adjusts forward and reverse flow to clean the endoscope lens when you tap the footswitch
- Six cycle settings allow you to customize performance
- Auto-prime prohibits drips or air bubbles from forming over the lens
- Easier set-up and switching between endoscopes with pump head design
- Covers a broad range of endoscopes by multiple manufacturers

*Fogging and accumulated debris can decrease visualization and necessitate removing the scope during a surgical procedure for cleaning.*



## The Endo-Scrub® 2 and IPC® Difference

Significantly Improves the Ability to Operate in the Presence of Bleeding<sup>1-6</sup>



*Software-driven Endo-Scrub® 2 Lens Cleaning Sheaths ensure that fluid is evacuated from the sheath to keep the lens from fogging.*

## IPC® and Endo-Scrub® 2 System

- Kit\* includes:
  - 1898001 IPC® Console
  - 1897821 Power Cord, 6 Meter, IEC 320, 115V
  - 1898430 IPC Multi-Function Footpedal
  - 1898851 IPC Manual
  - 1897510 Basket
  - 1852000 Endo-Scrub® 2 Footswitch
  - 1991015 Endo-Scrub 2 Fingerswitch

IPCES2SYSKIT\* IPC and Endo-Scrub® 2 System

\*This kit is only available in the US. If you are not in the US, please order the above items individually.