




## Shepard Intrauterine Insemination Sets

Used for the introduction of washed spermatozoa into the uterine cavity. The rigidity of the inner cannula facilitates introduction of the catheter through the internal cervical os and into the uterus.

- The soft, flexible distal tip eases placement and promotes patient comfort.
- The more rigid proximal segment provides enhanced catheter control and positioning.
- The adjustable positioner allows optimum placement of the catheter within the uterus.



Global Product Number	Order Number	Fr	Length cm	Unit of Sale	Description
G16454	J-IUIE-541009	5.4	10	box of 20	curved tip; with inner stainless steel cannula
G16465	J-IUIE-542009	5.4	20	box of 20	curved tip; with inner stainless steel cannula
G16463	J-IUIE-542001	5.4	20	box of 20	curved tip; with inner stainless steel cannula

### Modified Shepard Intrauterine Insemination Catheters

G16469	J-IUIE-542031	5.4	20	box of 20	curved tip; without inner cannula
G16468	J-IUIE-542030	5.4	20	box of 20	straight tip; without inner cannula