


NeoGen AR Nail Instrument Set



Code	Product Description	Qty
269900	NeoGen AR Nail Instruments Set	
269910	NeoGen AR Nail Instruments Set (Empty)-PPSU	1

NeoGen AR Nail Instruments Tray 1



Code	Product Description	Qty
267190	Outer Sleeve for Distal Locking Screw	2
267200	Drill Sleeve 4.0	2
267210	Trocar, Φ 4.0mm	1
267220	Depth Gauge for Distal Locking Screw	1
267250	Guide Wire with Hook for Cap	1
267340	Drill Bit, Φ 4.0mm	2
267390	Tip Threaded Guide Wire, Φ 3.2 mm	2
267410	Combination Wrench	1


267530	Cannulated Awl	1
267540	Connecting Rod	1
267560	Guide Wire, Φ 3.2 mm	2
269150	Drill Bit, Φ 5.2mm	1
269160	Trochar Sleeve 8.1	1
269170	Drill Sleeve 5.2	1
269180	Trochar 8.1	1
269190	Ball Tip Guide Wire 2.5	1
269210	Wrench, SW5 for Bolt	1
269360	Drill Bit, Φ 4.3mm	2
269370	Drill Sleeve 4.3	2
269380	Wrench for Stopper	1
269390	Stopper, Φ 4.3mm	1

NeoGen AR Nail Instruments Tray 2




Code	Product Description	Qty
267120	Insertion Connecting Rod	1
267130	Insertion Connector	1
267230	Screwdriver for Distal Locking Screw	1
267240	Screwdriver, Hexagonal with Spherical Head	1
267260	Screwdriver for Cap, Cannulated	1
267270	Extended Rod, Cannulated	1
267290	Cleaning Stylet	1
267300	Extraction Screw for Nail	1
267320	Extraction Screw for AR Blade	1
267480	Key for AR Blade	1
267510	Slotted Hammer Guide for AR Blade	1
269140	Position Rod Drill, Φ 5.2mm	1
269200	Reducer	1
269350	Position Rod	1
269470	Guide Sleeve 17	1

NeoGen AR Nail Instruments Tray 3



Code	Product Description	Qty
267180	Direct Measuring Device for Guide Wire Φ 3.2mm	1
267330	Screwdriver, Hexagonal	1
267360	Reamer for AR Blade	1
267370	Cannulated Drill Bit for AR Blade	1
267400	Stopper, Φ 10mm	1
267470	Tissue Protector	1
267520	Radiographic Ruler	1
269240	Guide Wire Holder	1
269280	Reamer Shaft	2
269340	E Block	1
269430	T- Handle with Quick Coupling	1
269450	Entry Reamer, Φ 17.0mm	1
269460	Protection Sleeve	1

NeoGen AR Nail Instruments Tray 4



Code	Product Description	Qty
267140	Bolt for Aiming Arm 130°	1
267150	Protection Sleeve for AR Blade	1
267160	Trocar Sleeve 3.2	1
267170	Trocar, Φ 3.2mm	1
267310	Slotted Hammer	1
267350	Screwdriver for AR Blade	1
267450	Hammer	1
267460	Buttress Nut	1
269100	Proximal Aiming Bar for Long Nail	1
269110	Distal Aiming Bar for Long Nail	1
269120	Targeter	1
269250	Proximal Bolt	1
269260	Bolt for Targeter	1
269270	Distal Bolt	1
269290	Reamer, Φ 9mm	1
269300	Reamer, Φ 10mm	1
269310	Reamer, Φ 11mm	1
269320	Reamer, Φ 12mm	1
269330	Reamer, Φ 13mm	1
269400	Insertion Handle	1
269410	Aiming Arm 130° for Static Locking	1
269420	Connecting Screw	1
269440	Aiming Arm 130° for Dynamic Locking	1