



Indigo™ High-Speed Otologic Drill

VERSATILE. CONVENIENT. POWERFUL.



The new Indigo™ High-Speed Otologic Drill with straight and angled attachments offers excellent surgical access and visibility. Its small, lightweight body provides power and performance.

Indigo™ High-Speed Otologic Drill

VERSATILE. CONVENIENT. POWERFUL.

The new Indigo High-Speed Otologic Drill offers excellent surgical access and visibility in a wide range of otologic procedures.

Indigo™ Components

PRODUCT	DESCRIPTION	QTY
1845000	Indigo High-Speed Otologic Drill	1
1845010	Indigo Otologic Straight Attachment	1
1845020	Indigo Otologic Angled Attachment	1
1845040	Indigo Instrument Tray	1
1845030	Indigo Tubing	5



1845000



1845010



1845020



1845030

Indigo Drill System Key Features:

- ① High speed and ample torque for fast debulking of the mastoid
- ② Available with straight and angled attachments
- ③ Improved surgical access and visibility in narrow otologic spaces
- ④ Works with a wide variety of Medtronic Xomed otology burs
- ⑤ Grooves on attachment nose design allows for more stable irrigation
- ⑥ Powered by the IPC™ System with easy touchscreen menus



For further information, please contact your Medtronic ENT representative. You may also consult our website at www.MedtronicENT.com.

Medtronic ENT

6743 Southpoint Drive N
Jacksonville, FL 32216
USA
www.MedtronicENT.com
Toll free: (800) 874-5797
Fax: (800) 678-3995

International Telephone Numbers

Adriatic Regional Office 385-1-488-1120
Australia 1800-668-670
Baltic Regional Office 37-1-67560226
Belgium 32-2456-09-09
Canada 1800-217-1617
China 86-21-50800998
Czech Republic 420-2-9657-9580
France 33-470-679-800
Germany 49-2159-8149-209
Greece 30-210-67-79-099
Hong Kong 852-2919-1312
Hungary 36-30-5052987
India 91-22-26836733
Israel 972-9-972-4400
Italy 39-02-24137-324
Japan 81-3-6430-2017
Korea 82-2-3404-3600
Lebanon 961-1-370-670
Luxembourg 32-2456-09-09
Netherlands 31-45-566-8800

Poland 48-22-4656900
Russian Federation 7-495-580-73-77
Singapore 65-6776-6255
South Africa 27-11-466-1820
Spain 34-91-625-05-40
Taiwan 886-2-2183-6000
UK 44-1923-205-166
USA 1-904-296-9600

